



# START


STents And Radiation Therapy

**SUMMARY PRESENTATION**

# STents And Radiation Therapy (START)

**Purpose:** To assess the safety and effectiveness of intracoronary beta radiation using a <sup>90</sup>Sr / <sup>90</sup>Y source train following successful coronary intervention in patients with in-stent restenosis.


## STUDY SCHEMA




**Design:** Prospective, multi-center (50 sites in N. America & Europe), randomized, placebo-controlled, triple-masked clinical trial.

## CLINICAL OUTCOME ANALYSIS

Mean Lesion Length 16.1 mm ± 7.4



## Significantly Reduced Major Adverse Cardiac Events 360 days




## 8 MONTH SAFETY RESULTS

PARAMETER	PLACEBO	Sr-90
<b>Death</b>	1 (0.4%)	3 (1.2%)
<b>MI</b>	7 (3.0%)	4 (1.6%)
Q-wave	0	0
non-Q-wave	7	4
<b>Aneurysm<sup>1</sup></b>	0 (0%)	1 (0.5%)
<b>Thrombosis</b>	1 (0.4%)	0
In-hospital - 30 days	1	0
31 - 240 days	0	0 <sup>2</sup>
<b>Angiographic Total Occlusions</b>	7 (3%)	8 (3.3%)
Patients with new stent <sup>3</sup>	4/35	3/42
Patients with no new stent <sup>4</sup>	3/153	5/156


1. No new aneurysm formation: 1 patient with aneurysm present at baseline showed no significant change at follow-up
2. One patient recently adjudicated by CEC had thrombosis at day 244
3. 74% of patients received ≤ 60 days of adjunctive anti-platelet therapy
4. 88% of patients received ≤ 60 days of adjunctive anti-platelet therapy

# ANGIOGRAPHIC QCA ANALYSIS


## Segment Analysis



## 8 Month Late Loss




## Methodology



## LONG LESION SUBGROUP

Lesions  $\geq 15$  mm  
(Mean Lesion Length: 21.8 mm  $\pm$  5.3)



## Clinical Summary

- Pre-specified hypotheses were achieved with statistical significance
- TVF .....reduced by 31% (p=0.039)
  - MACE .....reduced by 31% (p=0.039)
  - TVR .....reduced by 34% (p=0.028)
  - TLR .....reduced by 42% (p=0.009)

## Angiographic Summary

- Pre-specified restenosis hypotheses were achieved with statistical significance
- Stent Segment .....reduced by 66% (p<0.001)
  - Analysis segment .....reduced by 36% (p=0.001)

## Safety Summary

- Placebo vs Sr-90
- No difference in Death or MI .....(8 vs 7)
  - No difference in Late Thrombosis .....(0 vs 1)
  - No difference in Total Occlusions .....(7 vs 8)
  - No difference in New Aneurysm Formation .....(0 vs 0)

## Long Lesion SubGroup Summary

Sr-90 significantly reduced restenosis rate and TVR in patients presented with long lesions.

## Conclusion

The Beta-Cath™ System has been shown to be safe and effective for the treatment of in-stent restenosis.

## STUDY ENDPOINTS

Primary Efficacy Endpoint:

- 8 month Target Vessel Failure (TVF)

Secondary Efficacy Endpoints:

- 8 month angiographic restenosis, in-stent MLD, and late loss

Safety Endpoints:

- 8 month MACE and aneurysm formation

## DOSIMETRY METHODS

- Reference vessel diameter (RVD) was determined visually after completion of coronary intervention
- Dose prescription point calculated @ 2 mm from center of source axis:

18.4\* Gy in RVD  $\geq$  2.7 -  $\leq$  3.3 mm

23\* Gy in RVD  $>$  3.3 -  $\leq$  4.0 mm

\*NIST dose, March 2000

- <sup>90</sup>Sr / <sup>90</sup>Y has a 28.8 year half-life and a short treatment time of 3 to 5 minutes.

## INCLUSION / EXCLUSION CRITERIA

Major Inclusion Criteria:

- Patients over 18 years of age
- Single lesion in single native coronary vessel (diameter 2.7 - 4.0 mm)
- In-stent restenosis  $>$  50% (by visual estimate)
- Lesion length treatable with 20 mm balloon

Major Exclusion Criteria:

- Multivessel coronary intervention
- Unsuccessful treatment ( $>$ 30% residual stenosis) of target lesion
- Recent ( $<$ 72 hours) MI
- LVEF  $<$  30%
- Unprotected left main disease
- Anticipated use of ReoPro® (Eli Lilly & Company) or placement of a second stent
- Prior chest radiotherapy

### PROCEDURE DETAILS

	Placebo	Sr-90
<b>Debulking Devices (%)</b>		
DCA	0.9	0.0
RA	39.8	43.9
ELCA	7.4	5.7
<b>New Stents* (%)</b>	<b>19.8</b>	<b>20.9</b>

\* "Bail-out" stent use was reserved for severe residual stenoses after radiation delivery.

### DEVICE PERFORMANCE

	Patients	Percent
Total Patients Enrolled	476	100.0%
Successful Treatment	467	98.1%
Catheter not cross lesion	6	1.3%
Sources not sent	3	0.6%

### BASELINE FINDINGS

	Placebo (n=232)	Sr-90 (n=244)
<b>Clinical Characteristics</b>		
Age (yrs)	61.1	61.5
Men (%)	63.4	68.4
Diabetes (%)	32.3	30.7
Smoking (%)	8.1	12.5
Prior MI (%)	47.8	46.7
Prior CABG (%)	23.7	21.4
<b>Angiographic Characteristics</b>		
Vessel Diameter, mm	2.77	2.76
MLD, mm	0.98	0.98
% Stenosis	64.2	64.2
Lesion Length, mm	16.0	16.3
% LAD	41.3	43.2

## Definitions

Target Vessel Failure (TVF) = Death attributed to the target vessel, MI and TVR

MACE = Death, MI, emergent CABG and TVR

Target Vessel Revascularization (TVR) = Any clinically-driven repeat percutaneous intervention of the target vessel or bypass surgery of the target vessel

Target Vessel Revascularization (TLR) = Any clinically-driven repeat percutaneous intervention of the target lesion or bypass surgery of the target vessel

**Best Vascular, Inc.**  
 4350 International Blvd, Suite A  
 Norcross, Georgia 30093  
 Tel: +1 770-717-0904  
 Fax: +1 770-717-1283

 Novoste, Beta-Cath and the Beta-Cath System design logo are trademarks of Best Vascular, Inc. U.S. Patent Numbers 5,683,345; 5,899,882 and 6,013,020. Other patents pending. © 2000 Best Vascular, Inc. All rights reserved

**International Use Only**  
**1-800-NOVOSTE**  
 (1-800-668-6783)  
[www.novoste.com](http://www.novoste.com)

MT 368-00 02/01